



SAFE-Q for iPad Operation Manual

December 23, 2019

1st Edition



The Japanese Society for Surgery of the Foot

Uninstalling the App	3
Launching the App for the First Time	4
Home	5
Registering Patient Information	6
Answering Questions	7
Completion-Checking Results	8
Checking Results	9
Comparing with Past Data	10
Mean Score + SD	11
Searching Patient Information	12
Deleting Patient Information	13
Editing Patient Information	14
Viewing Data List	15
Menu Display	16
Viewing Data on PC	17
Importing Registered Data to iPad	23


Uninstalling the App

1. Access the App Store and tap "Search" at the top right of the screen.
2. Enter "SAFE-Q" to search it.
3. Check the contents and tap "Get."
When "Get" changes to "FREE," tap the button another time.
*Enter your Apple ID and password.
4. Installation starts.

iPad is a registered trademark of Apple Inc. in the United States and other countries.

App Store is the name of the service that allows you to download application software for iPhone, iPad, iPod, and iPod touch operated by Apple Inc.


◆ iPad Home screen



Tap the "App Store" icon.



◆ App Store



Search for "SAFE-Q."

◆ Download app



Tap the **FREE** button.
When it changes to "FREE," tap the button again.

Launching the App for the First Time


After completing installation, tap the "SAFE-Q" icon to launch the app.

The Disclaimer / Terms of Service window will open.
After confirming the contents, tap "Agree and Use."

Go to the Home page.


* The "Disclaimer / Terms of Use" window is displayed only when the app is started for the first time, and will not be displayed from the second and subsequent launches.

◆ iPad Home screen




Tap the icon.

◆ Disclaimer / Terms of Use window



Tap "Agree and Use."

◆ Home page



Go to the Home page.
From the second time this app is started and onwards, this window will always be displayed first.


In this app, patients tap windows to answer the questionnaire themselves. Questionnaire results are calculated automatically. Each patient can take this questionnaire up to nine times.

Register the information of new patients from "Register." From the second time onwards, tap "Search patient" to open the patient information page.

To register new patients.
(See P6.)

◆ Home page

To check information on
registered patients. (See P12.)



Go to Arthrex.

Registering Patient Information

For new patients, enter at ID, Name, Furigana, Date of birth, Disease name, and Sex.
ID, Name, Furigana, and Date of birth are required fields, so be sure to enter them. Disease name and Sex are optional.

To go to questions, tap "Register and go to question" at the Check registered information window.
To register only, tap "Register." This ends registration.
To register another information, tap "Home" at the Registration completed window.

◆ Register window

Self-Administered Foot Evaluation Questionnaire (SAFE-Q)

Register patient information

ID : test01 (Required)
(Within 12 one-byte characters)

Name : test (Required)
(Within 24 one-byte characters)

Furigana : Enter furigana.

Date of birth : 01/01/1959 (Required)

Sex : Male Female

Disease name : Acquired deformation : Hallux valgus
Foot inflammation and foot disorders with systemic disease : Rheumatoid arthritis

Go to Confirm window

Home Register Search patient List

Enter patient information

After entering the required information, tap "Go to Confirm window."

* ID (Within 12 one-byte characters)
Name (Within 12 two-byte characters)

* To check if the ID is registered, After entering only "ID," tap "Go to Confirm window." If an already registered ID is entered, the corresponding registered patient information will be displayed.

◆ Check registered information window

Self-Administered Foot Evaluation Questionnaire (SAFE-Q)

Check registered information.

ID : test01

Name : test

Furigana :

Date of birth : 01/01/1959

Sex : Male

Disease name : Acquired deformation : Hallux valgus
Foot inflammation and foot disorders with systemic disease : Rheumatoid arthritis

Correct registered information Register Register and go to questions

Home Register Search patient List

Check registered information

"Correct registered information"
→Correct the entered information.

"Register"
→Register patient information only.

"Register and go to questions"
→Starts questions after registering.

◆ Registration completed window

Self-Administered Foot Evaluation Questionnaire (SAFE-Q)

Registration completed.

Go Home

Home Register Search patient List

◆ Start questions window (See P7.)

Self-Administered Foot Evaluation Questionnaire (SAFE-Q)

About This Survey

The following questionnaire will ask you about your feet, and the difficulties you may be experiencing or have in your daily life. This includes emotional difficulties in addition to pain and physical functions. This questionnaire will also ask you questions about quality of life (QOL), which may be affected by foot diseases and injuries.

Your candid feedback will contribute to future treatments and support for feet problems.

Thank you for your cooperation in this questionnaire.

Foot here refers to the area enclosed by the square in the figure on the right, from the shin to the toes, and does not include the knee.

== Note on your replies ==

① For each question, please try to remember your condition this last week or this last month.

② The instructions for answering the questions are explained for each question. Please read them carefully before replying. They will take about 10 minutes although there are individual differences.

③ In this questionnaire, there are two ways of answering questions as follows.

I. Write ✓ in the □ of your answer. OR

II. Write x on the line.

Start of Questions

Home Register Search patient List

Answering Questions (by Patients)

Tap "Start questions" to start the questions (Please hand the iPad to the patient)

After selecting an answer, tap "Next" to go to the next question.

* Tap "Previous" to return to the previous question and answer again.

Tap the "Japanese Society for Surgery of the Foot." logo at the top right of the window to return to the Home page.

"Note" If you return to the Home page in the middle of a question, your answer will not be saved, and you will need to answer the question again.

There are two ways of answering questions.

I. Writing ✓ in the ☐ of your answer. OR

II. Write x on the line.

◆ Start questions window

Self-Administered Foot Evaluation Questionnaire (SAFE-Q)

About This Survey

The following questionnaire will ask you about your feet, and the difficulties you may be experiencing or have in your daily life. This includes emotional difficulties in addition to pain and physical functions. This questionnaire will also ask you questions about quality of life (QOL), which may be affected by foot diseases and injuries.

Your candid feedback will contribute to future treatments and support for feet problems.

Thank you for your cooperation in this questionnaire.

"Foot" here refers to the area enclosed by the square in the figure on the right, from the shin to the toes, and does not include the knee.

== Note on your replies ==

① For each question, please try to remember your condition during the past month.

② The instructions for answering the questions are explained below. Please read them carefully before replying. They will take about 10 minutes although the actual answering time is shorter.

③ In this questionnaire, there are two ways of answering questions.

I. Write ✓ in the ☐ of your answer. OR

II. Write x on the line.

Start of Questionnaire ▶

◆ Questions window (Answering by ticking ✓)

Self-Administered Foot Evaluation Questionnaire (SAFE-Q)

Q2

Have you had difficulty in sleeping due to foot pain in the past week?

(select the one that applies and place an ✓ in the corresponding box)

Never ☐ Occasionally ☐ Sometimes ☒ Often ☐ Always ☐

◀ Back Next ▶

Tapping "Previous" returns to the previous question.

Tapping "Next" proceeds to the next question.

Tap the logo to return to the Home page.

◆ Questions window ((Answering by adding a cross on the line x))

Self-Administered Foot Evaluation Questionnaire (SAFE-Q)

問3

How intense was the most severe pain you experienced in your feet in the past week?

(select the one that applies and place an "x" in the corresponding box)

0 ————— 10

No pain Worst pain imaginable

Touch on the line.

◀ Back Next ▶


Touch the line to show the degree of pain.

Completion-Checking Results


When the patient completes all questions, have the patient return the iPad to you.
Tap the “Japanese Society for Surgery of the Foot” logo at the bottom right of the window to go to the Results window.

* To prevent patients from looking at their own results, the “See results” button will not be displayed

◆End questions End questions window




◆Results window



Checking Results

The latest two answers registered are displayed. You can select the answers of any questionnaires registered and compare them. Just select the registered questionnaires (registration number) you want to compare and tap “Show results.” (e.g.: The figure below compares the values of the second and third questionnaires registered.)

◆ Results window



Select the registered answers you want to compare and tap “Show results.”


To clear the window, tap “Clear.”

Tap “Compare with past data (Radar chart)” to go to the radar chart window.

Comparing with Past Data

The latest two answers registered are displayed. You can select the answers of any questionnaires registered and compare them. Just select the registered questionnaires (registration number) you want to compare and tap "Show results." The radar chart is displayed. (e.g.: The figure below compares the values of the second and third questionnaires registered.)

◆Results window (Radar chart)



To change the registered answers you want to compare, select the number of the registered answer and tap "Show results."

To clear the window, tap "Clear."


Tap "Results list" to return to the Results window.

Tap "Mean score + SD" to go to the mean score +standard deviation chart window.

Mean Score + SD

This is the mean score + SD of the latest answers of all patients registered on the iPad.
(e.g.: The figure below shows the mean score + SD of the answers from patient A's third questionnaire, answers from patient B's second questionnaire, and answers from patient C's first questionnaire.)

◆ Mean score + SD



Tap "Search patient" to return to the Search patient window. (See P12.)


Searching Patient Information

Patient information is listed on the Search patient window in the order registered.

When searching for a specific patient, select the item to use for the search from "ID," "Name," "Furigana," "Date of birth," and "Disease name." Enter a keyword and tap the "Search" button to display the search results.

From this window, you can go to the "Question," "Patient Information," and "Show results" windows of the patient. You can also delete information.

◆ Search patient window



Select the item to use for the search, enter a keyword, and tap the "Search" button.

To search again, tap the "Clear" button.


- Question (See P7.)
Goes to the Question window.
- Patient information (See P6.)
Goes to the Check patient information window.
- Show results (See P9.)
Goes to the Show results window.
- Delete (See P13.)
Deletes the registered patient information.

Tap "Delete all" to delete the information of all registered patients. (See P13.)

Deleting Individual or All Patient Information


To delete individual patient information, tap "Delete." The confirmation window will appear twice. Tap "Yes" on each window to delete the information. Select "No" if you do not want to delete the information.

[Note] Once deleted, information cannot be restored. For deleted information, you need to re-register the patient and the patient needs to answer the questionnaire again.



To delete all patient information, tap "Delete all." The confirmation window will appear twice. Tap "Yes" on each window to delete all the information. Select "No" if you do not want to delete the information.

[Note] Once deleted, information cannot be restored. For deleted information, you need to re-register the patient and the patient needs to answer the questionnaire again.




Editing Patient Information

Tap the “Patient Information” of the patient you want to edit from the patient search list. The Confirm registered information window appears.

From this window, you can go to the windows for deleting data, editing registered information, and the Questions window.

[Note] Once deleted, information cannot be restored. For deleted information, you need to re-register the patient and the patient needs to answer the questionnaire again.

Tap “Edit data” to go to the Edit registration window.
Tap the item to be changed, re-enter the information, and then tap the “Change” button to proceed to the Completed window.



◆ Check registered information window

Self-Administered Foot Evaluation Questionnaire (SAFE-Q)

Registered information.

ID : test01

Name : test

Furigana :

Date of birth : 01/01/1959

Sex : Male

Disease name : Hallux valgus
Rheumatoid arthritis

Date of register : 1st 11/23/2019 / 2nd 11/23/2019 / 3rd 11/25/2019

Delete data Edit data Go to questions (4th time)

Home Register Search patient List

Tap “Edit data.”

◆ Edit registered information window

Self-Administered Foot Evaluation Questionnaire (SAFE-Q)

Re-enter the items to be changed.

ID test01 (Required)

(Within 12 one-byte characters)

Name test (Required)

(Within 24 one-byte characters)

Furigana Enter furigana.

Date of birth 12/31/1959 (Required)

Sex Male Female

Disease name Acquired deformation : Hallux valgus
Foot inflammation and foot disorders with systemic disease :
Rheumatoid arthritis

Enter disease name.

Cancel Change

Home Register Search patient List

Tap the information to be changed, re-enter, and tap “Change” to end.

◆ Completed window

Self-Administered Foot Evaluation Questionnaire (SAFE-Q)

Changes completed

ID : test01

Name : test

Furigana :

Date of birth : 12/31/1959

Sex : Male

Disease name : Acquired deformation : Hallux valgus
Foot inflammation and foot disorders with system...

Return to Search patient window

Home Register Search patient List

Viewing Data List


On the Data list window, you can check the date of registration of the patient information and basic information such as the patient ID, and the date the questionnaire answers were registered.

The patient list can be sorted in ascending or descending order by "Registration date," "ID," "Date of birth," and "Sex."

Tap "List" to list the registration dates of answers and answer data.

Select the item you want to sort and the descending or ascending order. Tap the "Sort" button.


◆ Data list window



Tapping "List" changes to the Individual answer list.

Tap the "Return" button to return to the Data list window.

◆ List window (List of answers)



To check answers, swipe the List window to the left.


Menu Display

From the menu at the bottom of the window, you can go to the "Home," "Register," "Search patient," and "List" windows.
* As the Questions window are operated by the patients, there is no menu.


◆ Home page





◆ Search patient window



◆ Register window



◆ Data list



Viewing Data on PC

The contents of the data list can be output in CSV format which can be opened on PCs. (Download with iTunes)

1. Tap "Create output file" at the Data list window.
2. Select whether to set a password the output file.
3. Connect the iPad to the PC and start iTunes.
4. Selecting "App" in iTunes shows the SAFE-Q app icon in "File Sharing."
5. CSV data whose name contains the date it was created will be displayed. Download it to the desktop by dragging and dropping. Or, specify the folder to save it where you want.

[Note] CSV data for which passwords can be set are files for writing and are not master data.

◆Data list window

Self-Administered Foot Evaluation Questionnaire (SAFE-Q)

Data list

Date of registration	ID	Name	No. of patients
11/23/2019	test01	test	12.
11/25/2019	test02	test02	02.
			03.

Enter the password and save the output file to lock the file with the password.

Connect to computer with USB cable.

Self-Administered Foot Evaluation Questionnaire (SAFE-Q)

Data list

Please enter your password.

OK

Save without setting password.

Cancel


Please enter your password.

OK

Save without setting password.


Cancel

◆iTunes screen



Viewing Data on PC 2

◆ iTunes window



Download the file whose name contains the date to the desktop by drag and drop.

Or specify the folder to save it.

"PatientData_SAFE-Q.csv" is the master data. (Automatically created file.)

Viewing Data on PC 3: Downloaded CSV File

Open the downloaded CSV file. (Can be opened using Excel.)
You can check the data of the patients registered in iPad (Date of registration, ID, Name, Furigana, Date of birth, Sex, Disease name), date of registration of question, and the answers given to questions.

◆Downloaded CSV file (The following is an example of a window opened using Excel.)

The screenshot shows an Excel spreadsheet with the following data:

Paste patient information										Subscale	Pain and Pain-Related								Shoe-related
Patient name	Furigana	Date of birth	Gender	Disease name 1 (Large category)	Disease name 1 (Subcategory)	Disease name 2 (Large category)	Disease name 2 (Subcategory)	Disease name 3 (Large category)	Disease name 3 (Subcategory)	Question No.	Q1	Q2	Q3	Q4	Q5	Q6	Q7	Q8	Q9
test		12/31/1959	Male	Acquired deformation	Hallux valgus	Acquired deformation	Rheumatoid arthritis			11232019	2	2	6.7	2	2	1	1	3	3
test02		02/02/1940	Female	Congenital deformation	Congenital clubfoot	Acquired deformation	Flat feet						0						
test03		03/03/1933	Male	Trauma	Talar fracture	Paralytic foot	Flaccid paralysis						0						


Green callout boxes with labels:

- Patient data (points to the first four columns of the patient information section)
- Date of registration of question (points to the 'Question No.' column)
- Answers to questions (points to the 'Q1' through 'Q9' columns)

Viewing Data on PC 3: How to Use CSV Files

Answers to questions can be used by cutting and pasting to the "[SAFE-Q App]_data aggregation" excel files.

◆[SAFE-Q App]_Data aggregation.xls



Questionnaire results pasting

Graph data

Graph

Made up of three sheets.

[Questionnaire results pasting] sheet

For pasting the data of the output CSV file (patient information and answers).

[Graph data] sheet

After pasting the data on the [Questionnaire results pasting] sheet, you can check the data aggregated for the graph on this sheet.

You can check the subscale total, total score, average, number of people (e.g.), standard deviation, standard error, and search for disease names.

[Graph] sheet

Mean score + SD graph can be checked.

Viewing Data on PC 5: Using Excel for Aggregating [SAFE-Q App]_Data

Copy the contents of the downloaded CSV file and paste to [SAFE-Q app] _Data aggregation excel file.

After selecting the target data in the CSV file, copy ([Right click → Copy]), paste the patient information is in the column for pasting patient information (green column in the figure below), and the answers in the column for pasting answers (blue in the figure below) ([Right click → Paste]).

◆Downloaded CSV file

Downloaded CSV file

	C	D	E	F	G	H	I	J	K	L	M	
1	Patient na	Furigana	Date of bir	Gender	Disease n	Disease n	Disease n	Disease n	Disease n	Disease n	Question I	Q1
2	test		12/31/195	Male	Acquired c	Hallux valg	Foot inflar	Rheumatoid arthritis			11232019	
3	test02		02/02/194	Female	Congenita	Congenita	Acquired c	Flat feet				
			03/03/193	Male	Trauma	Talar fract	Paralytic f	Flaccid paralysis				

Copy & Paste

	M	N	O	P	Q	R	S	T	U	V	W	X	
ase n	Disease n	Question I	Q1	Q2	Q3	Q4	Q5	Q6	Q7	Q8	Q9	Q10	Q11
thriti		11232019		2	2	6.7	2	2	1	1	3	3	2
						0							
ysis						0							

Copy & Paste

◆[SAFE-Q app]_Data aggregation Excel file Sheet name[For pasting questionnaire results]

Paste patient information																			
Patient name	Furigana	Date of birth	Gender	Disease name 1 (Large category)	Disease name 1 (Subcategory)	Disease name 2 (Large category)	Disease name 2 (Subcategory)	Disease name 3 (Large category)	Disease name 3 (Subcategory)	Subscale	Pain and Pain-related								Shoe-related
										Question No.	Q1	Q2	Q3	Q4	Q5	Q6	Q7	Q8	Q9
test		12/31/1959	Male	Acquired deformation	Hallux valgus	Foot disorders with	Rheumatoid arthritis			11232019	2	2	6.7	2	2	1	1	3	
test02		02/02/1940	Female	Congenital deformation	Congenital clubfoot	Acquired deformation	Flat feet						0						
test03		03/03/1933	Male	Trauma	Talar fracture	Paralytic foot	Flaccid paralysis						0						

Viewing Data on PC 6: Searching Patient Name

To check specific patient data only, search and extract the data.

◆[SAFE-Q app]_Data aggregation.xls Sheet name: [Graph data]

Paste patient information				
Patient name	Furigana	Date of birth	Gender	Disease name 1 (Large category)
test		12/31/1959	Male	Acquired deformation
test02		02/02/1940	Female	Congenital deformation
test03		03/03/1933	Male	Trauma

①Open the drop list of the disease name row.

Paste patient information

Date of birth Gender Disease name 1 (Large category)

12/31/1959 Sort A to Z Sort Z to A

02/02/1940 Sort by Color

03/03/1933 Clear Filter From "Disease name 1 (Large category)"

Filter by Color

Text Filters

Search

☒ (Select All)

☒ Acquired deformation

☒ Congenital deformation

☒ Trauma

☒ (Blanks)

OK Cancel

Paste patient information						
Patient name	Furigana	Date of birth	Gender	Disease name 1 (Large category)	Disease name 1 (Subcategory)	Disease name 1 (Large category)
test		12/31/1959	Male	Acquired deformation	Hallux valgus	nd foot

③The patient corresponding to the search results is displayed.
*The above figure shows the patient extracted when searched using "Rheumatism."

②Enter the disease name or copy-paste the disease name from the [Disease name list] and click OK.

Filter by Color

Text Filters

Search

☒ (Select All)

☒ Acquired deformation

☒ Congenital deformation

☒ Trauma

☒ (Blanks)

OK Cancel

To restore the file, check "Select all," and click the OK button to show all the data.

Importing Registered Data to iPad



You can share the “master data file (CSV data)” saved on your PC with your iPad.

1. Connect the iPad with the data to be transferred to your computer and start iTunes.
2. If you select “App” in iTunes, the SAFE-Q app icon will appear in “File Sharing.”
3. Drag and drop the saved CSV file to the “SAFE-Q document” field. Or, select the specified file from “Add File” and transfer it.

The message “Do you want to replace it with the item being moved?” appears. Select “Replace.”

[Note] Only master data can be transferred to the iPad. Since it will be overwritten, you should output the data once before transferring it (See P17).

◆ iTunes window




Drag the saved CSV file to the “SAFE-Q Document” area and drop it.

Or specify the file from “Add file.”

<Reference> How to Use the iPad


1. Tap

Tap the screen lightly with your finger. Do this to go to the next screen, etc.



2. Swipe

Slide your finger along the screen.




Move the x mark to the left and right while touching it to change its position.

* At this screen, you can also move the X mark by touching the line.


3. Drum roll (Pull down menu)

Move up and down to select the item.



4. Keyboard

Use to register patients and input free words.



Tapping the keyboard icon closes the keyboard.

Q1. The App Store does not show the app.

In the top right search window of the App Store screen, enter "SAFE-Q."
The app will be shown. (See P3.)

**Q3. I registered the wrong patient information.
What should I do?**

See page 14 of this manual and update the information. (See P14.)

**Q5. I deleted patient information by mistake.
What should I do?**

If you did not save the information in the output file, you will need to register the information again.
If you have saved the information in the output file, transfer the information to the app using iTunes. (See P23.)

* You cannot specify a patient in the saved output file and transfer the information of only that patient. Information can only be transferred all together.

Q2.App Store asks for a password when I press "Install." I don't know the password.

The password asked by App Store when installing the app is the password of your Apple ID.

**Q4. I forgot the password of the output file.
What should I do?**

Please output the file and set the password again. (See P17.)

For other questions on the app, please contact the email on the right.
(*For inquiries regarding the contents of the app, contact the Association.)

Dentsu, Sudler & Hennessey Inc.
2F Tsukiji Eto Bldg. 12-6 Tsukiji 1-chome, Chuo-ku, Tokyo 104-8427
For inquiries on the app: app-info@dsh.co.jp